

Product Datasheet

COX IV antibody (orb421032)

Description	Rabbit polyclonal antibody to COX IV
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Synthesized Peptide
Target	Cox4i1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	19530
Uniprot ID	P19783
Dilution Range	WB: 1:2000
Expiration Date	12 months from date of receipt.