

Product Datasheet

**beta Catenin CTNNB1 Antibody (monoclonal, 1F6)
(orb421109)**

Description

beta Catenin CTNNB1 Antibody (monoclonal, 1F6)

Species/Host

Mouse

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested Applications

FC, ICC, IF, IHC, WB

Immunogen

E. coli-derived human beta Catenin recombinant protein (Position: A2-K233).

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

Add 0.2ml of distilled water will yield a concentration of 500ug/ml.

Isotype

Mouse IgG1

Clonality

Monoclonal

MW

95 kDa

Uniprot ID


[P35222](#)

Dilution Range


Western blot, 0.1-0.5µg/ml Immunohistochemistry (Paraffin-embedded Section), 0.5-1µg/ml Immunocytochemistry/Immunofluorescence, 2µg/ml Flow Cytometry, 1-3µg/1x10⁶ cells

Expiration Date

12 months from date of receipt.




Flow Cytometry analysis of SiHa cells us...



Human Hela (Human COX2) Human HepG2 (Human Hsp70) Human SK-N-SH (Human GAPDH) Human L937 (Human IgG1)

WB analysis of beta Catenin using anti-b...



IF analysis of beta Catenin using anti-b...