



Product Datasheet

Elafin/Skalp PI3 Antibody (monoclonal, C7) (orb421110)

Description	Elafin/Skalp PI3 Antibody (monoclonal, C7)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	E. coli-derived human Elafin/Skalp recombinant protein (Position: A61-Q117).
Form/Appearance	Lyophilized
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
Storage	Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.
Note	For research use only
Application notes	Add 0.2ml of distilled water will yield a concentration of 500ug/ml.
Isotype	Mouse IgG1
Clonality	Monoclonal
Uniprot ID	P19957
Dilution Range	Immunohistochemistry (Paraffin-embedded Section), 0.5-1µg/ml ELISA (Cap), 1-5µg/ml
Expiration Date	12 months from date of receipt.



Sandwich
ELISA -
Recombinant
human
Elafi...



IHC analysis
of Elafin
using anti-
Elafin...