



Product Datasheet

Cardiac FABP Fabp3 Antibody (monoclonal, 6C4) (orb421121)


Description	Cardiac FABP Fabp3 Antibody (monoclonal, 6C4)
Species/Host	Mouse
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	E. coli-derived rat Cardiac FABP recombinant protein (Position: A2-A133).
Form/Appearance	Lyophilized
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
Storage	Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.
Note	For research use only
Application notes	Add 0.2ml of distilled water will yield a concentration of 500ug/ml.
Isotype	Mouse IgG1
Clonality	Monoclonal
MW	15 kDa
Uniprot ID	P07483
Dilution Range	Western blot, 0.1-0.5µg/ml Immunohistochemistry (Paraffin-embedded Section), 0.5-1µg/ml ELISA (Cap), 1-5µg/ml
Expiration Date	12 months from date of receipt.



WB analysis of Cardiac FABP using anti-C...



Sandwich ELISA - Recombinant rat Cardiac...



IHC analysis of Cardiac FABP using anti-...