



Product Datasheet

TPX2 antibody (orb421962)


Description	Mouse monoclonal antibody to TPX2
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IHC-P, IP, WB
Immunogen	Human TPX2
Target	TPX2
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Uniprot ID	Q9ULW0
Entrez	22974
Expiration Date	12 months from date of receipt.



Immunoprecipitation of TPX2 from nuclear...



Western blotting analysis of TPX2 using



Immunocytochemistry of TPX2 using monoclonal...