



Product Datasheet

CD3 antibody (orb421969)


Description	Mouse monoclonal antibody to CD3
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, ICC
Target	CD3
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 1-4 µg/ml.
Isotype	Mouse IgG2b
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Uniprot ID	P07766
Entrez	916
Dilution Range	Flow cytometry: Recommended dilution: 1-4 µg/ml.
Expiration Date	12 months from date of receipt.



Separation of human CD3 positive lymphoc...



Flow cytometry surface staining pattern ...



Immunocytochemistry staining of human CD...