

Product Datasheet

E-4031 (orb422286)

Description	E-4031 is a benzenesulfonamide antiarrhythmic agent; blocks the ATP-sensitive potassium channel.
Note	For research use only
Formula	$C_{21}H_{29}Cl_2N_3O_3S$
MW	474.4443
CAS Number	[113559-13-0]
SMILES	<chem>CS(=O)(NC1=CC=C(C(C2CCN(CCC3=NC(C)=CC=C3)CC2)=O)C=C1)=O.[H]Cl.[H]Cl</chem>
Expiration Date	12 months from date of receipt.