

Product Datasheet

AZD4635 (orb422434)

Description	AZD4635 is a novel adenosine 2A receptor (A2AR) inhibitor with a K_i of 1.7 nM.
Note	For research use only
Formula	$C_{15}H_{11}ClFN_5$
MW	315.7327
CAS Number	[1321514-06-0]
SMILES	<chem>C1C([H])=C(C([H])=C(C([H])([H])[H])N=1)C1=C(C2C([H])=C([H])C(=C([H])C=2[H])F)N=C(N([H])[H])N=N1</chem>
Expiration Date	12 months from date of receipt.