


Product Datasheet

CTGF antibody (orb422483)


Description	Goat polyclonal antibody to CTGF
Species/Host	Goat
Reactivity	Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC
Target	Connective Tissue Growth Factor
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Note	For research use only
Application notes	IHC: Paraffin embedded Rat Liver. Recommended concentration: 6-8µg/ml. Notes: Immunofluorescence: Strong expression of the protein seen in the cytoplasm of NIH3T3 cells. Recommended concentration: 5µg/ml.
Protein Sequence	SLYYRKMVGDMA
Clonality	Polyclonal
MW	38.1
Entrez	1490
Dilution Range	ELISA: 1:1000, IHC-P: 5 µg/ml
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of paraforma...



Immunohistochemical staining of Rat Live...



Immunohistochemical staining of Rat Live...