


Product Datasheet

F13B antibody (orb422592)

Description	F13B antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human Coagulation factor XIII B chain protein (260-403AA)
Target	F13B
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P05160
Dilution Range	WB:1:500-5000
Expiration Date	12 months from date of receipt.



Western blot analysis of 293 whole cell...