

Product Datasheet

SRRM4 antibody (orb422607)

Description	SRRM4 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant Human Serine/arginine repetitive matrix protein 4 protein (400-479AA)
Target	SRRM4
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	A7MD48
Dilution Range	IHC-P:1:500-1000
Expiration Date	12 months from date of receipt.