


Product Datasheet

KIR3DS1 antibody (orb422666)

Description	KIR3DS1 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF
Immunogen	Recombinant Human Killer cell immunoglobulin-like receptor 3DS1 protein (198-335AA)
Target	KIR3DS1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q14943
Dilution Range	IF:1:50-200
Expiration Date	12 months from date of receipt.



Immunofluorescence staining of HepG2 cel...