


Product Datasheet

ompA antibody (orb422682)

Description	ompA antibody
Species/Host	Rabbit
Reactivity	Bacteria
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Rickettsia conorii Outer membrane protein A protein (1734-2021AA)
Target	ompA
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q52657
Dilution Range	WB:1:500:5000
Expiration Date	12 months from date of receipt.



Western blot analysis of recombinant omp...