


Product Datasheet

CYP26C1 antibody (orb422693)

Description	CYP26C1 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF
Immunogen	Recombinant Human Cytochrome P450 26C1 protein (331-460AA)
Target	CYP26C1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q6V0L0
Dilution Range	IF/ICC:1:50:200
Expiration Date	12 months from date of receipt.



Immunofluorescence staining of HepG2 cel...