



## Product Datasheet

### TBC1D1 antibody (orb422708)

|                            |                                                                   |
|----------------------------|-------------------------------------------------------------------|
| <b>Description</b>         | TBC1D1 antibody                                                   |
| <b>Species/Host</b>        | Rabbit                                                            |
| <b>Reactivity</b>          | Human                                                             |
| <b>Conjugation</b>         | Unconjugated                                                      |
| <b>Tested Applications</b> | ELISA, IHC                                                        |
| <b>Immunogen</b>           | Recombinant Human TBC1 domain family member 1 protein (566-718AA) |
| <b>Target</b>              | TBC1D1                                                            |
| <b>Preservatives</b>       | Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol       |
| <b>Form/Appearance</b>     | Liquid                                                            |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.     |
| <b>Note</b>                | For research use only                                             |
| <b>Isotype</b>             | IgG                                                               |
| <b>Clonality</b>           | Polyclonal                                                        |
| <b>Uniprot ID</b>          | <a href="#">Q86TI0</a>                                            |
| <b>Dilution Range</b>      | IHC-P:1:500:1000                                                  |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                   |



Immunohistochemistry of paraffin-embedde...



Immunohistochemistry of paraffin-embedde...