



Product Datasheet

WDFY2 antibody (orb422736)

Description	WDFY2 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Immunogen	Recombinant Human WD repeat and FYVE domain-containing protein 2 protein (2-114AA)
Target	WDFY2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q96P53
Dilution Range	WB:1:500:5000, IF/ICC:1:50:200
Expiration Date	12 months from date of receipt.



Immunofluorescence staining of HeLa cell...



Western blot analysis of MCF-7 whole ce...