



Product Datasheet

CLCC1 antibody (orb422746)


Description	CLCC1 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, IP, WB
Immunogen	Recombinant Human Chloride channel CLIC-like protein 1 protein (446-551AA)
Target	CLCC1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q96S66
Dilution Range	WB:1:500:5000, IHC-P:1:200:500, IF/ICC:1:50:200, IP:1:200:2000
Expiration Date	12 months from date of receipt.



Immunohistochemistry of paraffin-embedde...



Western blot aanalysis of Colo320 whole ...



Immunoprecipitation analysis of HepG2 wh...