



Product Datasheet

INTS13 antibody (orb422754)

Description	INTS13 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IP, WB
Immunogen	Recombinant Human Integrator complex subunit 13 protein (573-706AA)
Target	INTS13
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9NVM9
Dilution Range	WB:1:500:5000, IP:1:200:2000
Expiration Date	12 months from date of receipt.



Immunoprecipitation
analysis of HeLa
who...



Western blot
analysis of HeLa
whole cel...