



Product Datasheet

NDUFAF3 antibody (orb422765)

Description	NDUFAF3 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	Recombinant Human NADH dehydrogenase [ubiquinone] 1 alpha subcomplex assembly factor 3 protein (91-184AA)
Target	NDUFAF3
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9BU61
Dilution Range	IHC-P:1:20:200
Expiration Date	12 months from date of receipt.



Immunohistochemistry of paraffin-embedde...



Immunohistochemistry of paraffin-embedde...