


Product Datasheet

SLC45A2 antibody (orb422792)

Description	SLC45A2 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human Membrane-associated transporter protein (1-46AA)
Target	SLC45A2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9UMX9
Dilution Range	WB:1:500:5000
Expiration Date	12 months from date of receipt.



Western blot analysis of MCF-7 whole cell lysate.