

Product Datasheet

MKL2 antibody (orb422796)

Description	MKL2 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant Human MKL/myocardin-like protein 2 protein (164-313AA)
Target	MKL2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9ULH7
Dilution Range	WB:1:500:5000, IHC-P:1:500:1000
Expiration Date	12 months from date of receipt.