


Product Datasheet

PLAGL2 antibody (orb422797)

Description	PLAGL2 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human Zinc finger protein PLAGL2 protein (324-469AA)
Target	PLAGL2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9UPG8
Dilution Range	WB:1:500:5000
Expiration Date	12 months from date of receipt.



Western blot analysis of HeLa whole cell...