


Product Datasheet

SPO11 antibody (orb422800)

Description	SPO11 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human Meiotic recombination protein SPO11 protein (28-109AA)
Target	SPO11
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9Y5K1
Dilution Range	WB:1:500:5000
Expiration Date	12 months from date of receipt.



Western blot analysis of Jurkat whole C...