

Product Datasheet

MAFB antibody (orb422801)

Description	MAFB antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human Transcription factor MafB protein (168-323AA)
Target	MAFB
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9Y5Q3
Dilution Range	WB:1:500:5000
Expiration Date	12 months from date of receipt.