



Product Datasheet

LILRB3 antibody (orb422826)

Description	LILRB3 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Immunogen	Recombinant Human Leukocyte immunoglobulin-like receptor subfamily B member 3 protein (470-624AA)
Target	LILRB3
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	075022
Dilution Range	WB:1:500:5000, IF/ICC:1:50:200
Expiration Date	12 months from date of receipt.



Immunofluorescence staining of NIH/3T3 C...



Western blot analysis of HepG2 whole ce...