



Product Datasheet

NAGPA antibody (orb422834)

Description	NAGPA antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant Human N-acetylglucosamine-1-phosphodiester alpha-N-acetylglucosaminidase protein (327-438AA)
Target	NAGPA
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9UK23
Dilution Range	WB:1:500:5000, IHC-P:1:500:1000
Expiration Date	12 months from date of receipt.



Immunohistochemistry of paraffin-embedde...



Western blot aanalysis of SH-SY5Y whole ...