



Product Datasheet

OPN1MW antibody (orb422839)

Description	OPN1MW antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Immunogen	Recombinant Human Medium-wave-sensitive opsin 1 protein (1-52AA)
Target	OPN1MW
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P04001
Dilution Range	WB:1:500:5000, IF/ICC:1:50:200
Expiration Date	12 months from date of receipt.



Immunofluorescence staining of HepG2 cel...



Western blot aanalysis of MCF-7 whole ce...