


Product Datasheet

POU1F1 antibody (orb422846)

Description	Rabbit polyclonal antibody to POU1F1
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P
Immunogen	Recombinant Human Pituitary-specific positive transcription factor 1 protein (34-130AA)
Target	POU1F1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P28069
Dilution Range	IHC-P:1:500:1000
Expiration Date	12 months from date of receipt.



Immunohistochemistry of paraffin-embedde...