


Product Datasheet

TBX4 antibody (orb422858)

Description	TBX4 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human T-box transcription factor TBX4 protein (355-511AA)
Target	TBX4
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P57082
Dilution Range	WB:1:500:5000
Expiration Date	12 months from date of receipt.



Western blot
analysis of
A549 whole
cel...