

Product Datasheet

PDE6C antibody (orb422871)

Description	PDE6C antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Immunogen	Recombinant Human Cone cGMP-specific 3',5'-cyclic phosphodiesterase subunit alpha protein (285-451AA)
Target	PDE6C
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P51160
Dilution Range	WB:1:500:5000, IF/ICC:1:50:200
Expiration Date	12 months from date of receipt.