


Product Datasheet

KIR2DS5 antibody (orb422888)

Description	KIR2DS5 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF
Immunogen	Recombinant Human Killer cell immunoglobulin-like receptor 2DS5 protein (162-242AA)
Target	KIR2DS5
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q14953
Dilution Range	IF/ICC:1:50:200
Expiration Date	12 months from date of receipt.



Immunofluorescence staining of HepG2 cel...