


Product Datasheet

POF1B antibody (orb422947)

Description	POF1B antibody
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human Protein POF1B protein (532-589AA)
Target	POF1B
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8WVV4
Dilution Range	WB:1:500:5000
Expiration Date	12 months from date of receipt.



Western blot analysis of K562 whole cell...