


Product Datasheet

TRPM3 antibody (orb422982)

Description	TRPM3 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF
Immunogen	Recombinant Human Transient receptor potential cation channel subfamily M member 3 protein (1606-1713AA)
Target	TRPM3
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9HCF6
Dilution Range	IF/ICC:1:50:200
Expiration Date	12 months from date of receipt.



Immunofluorescence staining of HeLa cell...