



Product Datasheet

Krt18 antibody (orb423027)


Description	Krt18 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant Mouse Keratin, type I cytoskeletal 18 protein (200-423AA)
Target	Krt18
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P05784
Dilution Range	WB:1:500:5000, IHC-P:1:20:200
Expiration Date	12 months from date of receipt.



Immunohistochemistry of paraffin-embedde...



Western blot aanalysis of Hela whole cel...



Immunohistochemistry of paraffin-embedde...