

Product Datasheet

RBP3 antibody (orb423041)

Description	RBP3 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Human Retinol-binding protein 3 protein (1146-1246AA)
Target	RBP3
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P10745
Dilution Range	WB:1:500:5000
Expiration Date	12 months from date of receipt.