


Product Datasheet

TH antibody (orb423071)

Description	TH antibody
Species/Host	Rabbit
Reactivity	Zebrafish
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Zebrafish Tyrosine hydroxylase protein (1-110AA)
Target	th
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9W7K7
Dilution Range	WB:1:500:5000
Expiration Date	12 months from date of receipt.



Western blot
analysis of
Positive WB
de...