


Product Datasheet

rbm24 antibody (orb423096)

Description	rbm24 antibody
Species/Host	Rabbit
Reactivity	Zebrafish
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Zebrafish RNA-binding protein 24 protein (2-78AA)
Target	rbm24
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q76LC6
Dilution Range	WB:1:500:5000
Expiration Date	12 months from date of receipt.



Western blot
analysis of
Positive WB
de...