


Product Datasheet

BHLHA15 antibody (orb423099)

Description	BHLHA15 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	Recombinant Human Class A basic helix-loop-helix protein 15 protein (2-72AA)
Target	BHLHA15
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q7RTS1
Dilution Range	IHC-P:1:20:200
Expiration Date	12 months from date of receipt.



Immunohistochemistry
of paraffin-embedde...