



Product Datasheet

HOOK3 antibody (orb423101)


Description	HOOK3 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, IP, WB
Immunogen	Recombinant Human Protein Hook homolog 3 protein (357-455AA)
Target	HOOK3
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q86VS8
Dilution Range	WB:1:500:5000, IHC-P:1:200:500, IP:1:200:2000
Expiration Date	12 months from date of receipt.



Immunohistochemistry of paraffin-embedde...



Western blot aanalysis of 293 whole cell...



Immunoprecipitation analysis of 293 whol...