

Product Datasheet

GALNT9 antibody (orb423102)

Description	GALNT9 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Human Adhesion G protein-coupled receptor G3 protein (98-244AA)
Target	ADGRG3
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q86Y34
Dilution Range	IHC-P:1:20:200
Expiration Date	12 months from date of receipt.