



Product Datasheet

APPL2 antibody (orb423106)

Description	APPL2 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC
Immunogen	Recombinant Human DCC-interacting protein 13-beta protein (399-512AA)
Target	APPL2
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q8NEU8
Dilution Range	IHC-P:1:500:1000
Expiration Date	12 months from date of receipt.



Immunofluorescence staining of HeLa cell...



Immunohistochemistry of paraffin-embedde...