


Product Datasheet

CLN8 antibody (orb423134)

Description	CLN8 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	Recombinant Human Protein CLN8 protein (247-286AA)
Target	CLN8
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9UBY8
Dilution Range	IHC-P:1:20:200
Expiration Date	12 months from date of receipt.



Immunohistochemistry
of paraffin-embedde...