

Product Datasheet

Mouse Anti-COX4I1(Mitochondrial Loading Control) antibody (orb423160)

Description

Mouse Anti-COX4I1(Mitochondrial

Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC-P, WB
Immunogen	Recombinant human COX4
Target	COX4I1(Mitochondrial Loading Control)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Monoclonal
MW	17 kDa
Uniprot ID	P13073
Dilution Range	WB=1:500-2000, IHC-P=1:100-500, IF=1:200-800, IHC-F=1:100-500, ICC=1:100-500
Expiration Date	12 months from date of receipt.