

## Product Datasheet

### CD3 antibody (orb423299)

|                          |                                               |
|--------------------------|-----------------------------------------------|
| <b>Description</b>       | Mouse monoclonal antibody to CD3              |
| <b>Species/Host</b>      | Rat                                           |
| <b>Reactivity</b>        | Mouse                                         |
| <b>Conjugation</b>       | Unconjugated                                  |
| <b>Note</b>              | For research use only                         |
| <b>Application notes</b> | Monoclonal Antibodies, Anti-Mouse Lymphocytes |
| <b>Clonality</b>         | Monoclonal                                    |
| <b>Expiration Date</b>   | 12 months from date of receipt.               |