


Product Datasheet

Human Fibrinogen antibody (orb4249)

Description	Human Fibrinogen antibody
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Fibrinogen from human plasma
Target	human Fibrinogen
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	240 kDa
CAS Number	9001-32-5
Uniprot ID	P02671
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.



Western blot: Lane 1: Human Plasma; Lane 2: ...