

Product Datasheet

Lassa GP1 Protein (orb426795)

Description	Recombinant of Lassa GP1 protein
Form/Appearance	Lassa GP1 protein solution contains 25mM Tris base and 25mM K ₂ CO ₃ in sterile filtered solution.
Storage	Lassa GP1 protein although stable at 4°C for 1 week, should be stored below -18°C. Please prevent freeze thaw cycles.
Tag	C-term 6x His
Note	For research use only
Application notes	Viral Antigens- Lassa Virus
Protein Sequence	205AA, partial.
Purity	Lassa GP1 is > 90% pure as determined by 10% PAGE (coomassie staining).
Source	Escherichia Coli
MW	23 kDa
Expiration Date	6 months from date of receipt.