

Product Datasheet

Human RXRA Protein (orb428240)

Description	RXRA Human Recombinant produced in E.Coli is a single, non-glycosylated polypeptide chain containing...
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Store at 4°C if entire vial will be used within 2-4 weeks. Store frozen at -20°C for longer periods of time. For long term storage, it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles.
Note	For research use only
Application notes	Recombinant & Natural Proteins
Protein Sequence	MLGLNGVLKV PAHPSGNMAS FTKHICAICG DRSSGKHYGV YSCEGCKGFF KRTVRKDLTY TCRDNKDCLI DKRQRNRCQY CRYQKCLAMG MKREAVQEER QRGKDRNENE VESTSSANE
Purity	Greater than 95.0% as determined by SDS-PAGE
MW	13.6 kDa
Expiration Date	6 months from date of receipt.