

Product Datasheet

Human SIGLEC5 Protein (orb428444)

Description

Recombinant of human SIGLEC5 protein

Preservatives

SIGLEC5 protein solution (0.25mg/ml) contains Phosphate Buffered Saline (pH 7.4).

Form/Appearance

Sterile Filtered colorless solution.

Storage

Stability: Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles

Note

For research use only

Application notes

Recombinant & Natural Proteins

Protein Sequence

ADLEKPVYEL QVQKSVTVQE GLCVLVPCSF
 SYPWRSWYSS PPLYVYWFRD GEIPYAEV
 ATNNPDRRVK PETQGRFRL GDVQKKNCSL
 SIGDARMEDT GSYFFRVERG RDVKYSYQQN
 KLNLEVTALI EKPDHFLEP LESGRPTRLS
 CSLPGSCEAG PPLTFSWTGN ALSPLDPETT
 RSELTTPR PEDHGTNLTC QMKRQGAQVT
 TERTVQLNVS YAPQTITIFR NGIALEILQN
 TSYLPVLEGQ ALRLLCDAPS NPPAHLWFQ
 GSPALNATPI SNTGILELRR VRSAEEGGFT
 CRAQHPLGFL QIFLNLSVYS LPQLLGPSCS
 WEAEGLHCRC SFRARPAPSL CWRLEEKPLE
 GNSSQGSFKV NSSSAGPWAN SSLILHGGLS
 SDLKVSCAW NIYGSQSGSV LLLQGRSNLG
 TGVVPAALLE PKSCDKTHTC PPCAPELLG
 GPSVFLFPPK PKDTLMISRT PEVTCVVVDV
 SHEDPEVKFN WYVDGVEVHN AKTKPREEQY
 NSTYRVVSVL TVLHQDWLNG KEYKCKVSNK
 ALPAIEKTI SKAKGQPREP QVYTLPPSRD
 ELTKNQVSLT CLVKGFYPSD IAVEWESNGQ
 PENNYKTPPP VLDSDGSEFFL YSKLTVDKSR
 WQQGNVFSCS VMHEALHNHY TQKSLSLSPG
 KHHHHHH

Purity

Greater than 90.0% as determined by SDS-PAGE.

Source

Sf9, Baculovirus cells

Expiration Date

12 months from date of receipt.