

Product Datasheet

Human SLAMF6 (Sf9) Protein (orb428472)

Description	Recombinant of human SLAMF6 (Sf9) protein
Preservatives	The SLAMF6 protein solution (0.5mg/ml) contains 1xPBS (pH 7.4) and 10% glycerol.
Form/Appearance	Sterile Filtered clear solution.
Storage	Stability: Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles
Note	For research use only
Application notes	Recombinant & Natural Proteins
Protein Sequence	ADPQSSLTPL MVNGILGESV TLPLEFPAGE KVNFIWLFN ETSLAFIVPH ETKSPEIHVT NPKQ GKRLNF TQSYSLQLSN LKMEDTGSYR AQISTKTS AK LSSYTLRILR QLRNIQVTNH SQLFQNMTCE LHLTCSVEDA DDNVSFRWEA LGNTLSSQPN LTVSWDPRIS SEQDYTCIAE NAVSNLSFSV SAQKLCEDVK IQYTDTKMHH HHHH
Purity	Greater than 95.0% as determined by SDS-PAGE.
Source	Baculovirus sf9 cells
Expiration Date	12 months from date of receipt.