

Product Datasheet

Human STAMPB Protein (orb428618)

Description	Recombinant of human STAMPB protein
Preservatives	STAMPB protein solution (1mg/ml) containing 20mM Tris-HCl buffer (pH 8.0), 0.15M NaCl, 10% glycerol and 1mM DTT.
Form/Appearance	Sterile Filtered colorless solution.
Storage	Stability: Store at 4°C if entire vial will be used within 2-4 weeks. Store, frozen at -20°C for longer periods of time. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Avoid multiple freeze-thaw cycles
Note	For research use only
Application notes	Recombinant & Natural Proteins
Protein Sequence	MGSSHHHHHH SSQLVPRGSH MGSMSDHGDV SLPPEDRVRA LSQLGSAVEV NEDIPRRRYF RSGVEIRMA SIYSEEGNIE HAFILYNKYI TLFIEKLPKH RDYKSAVIPE KKDTVKKLKE IAFPKAEELK AELLKRYTKE YTEYNEEKKK EAEELARNMA IQQELEKEKQ RVAQQKQQQL EQEQFHAFEE MIRNQELEKE RLKIVQEF GK VDPGLGGPLV PDLEKPSLDV FPTLTVSSIQ PSDCHTTVRP AKPPVVDRSL KPGALSNSSES IPTIDGLRHV VVPGR LCPQF LQLASANTAR GVETCGILCG KLMRNEFTIT HVLIPKQSAG SDYCNTEENEE ELFLIQDQQG LITLGWIH TH PTQT AFLSSV DLH THCSYQM MLPESVAIVC SPKFQETGFF KLTDHGLEEI SSCRQKGFHP HSKDPPLFCS CSHVTVDRA VTITDLR
Purity	Greater than 85% as determined by SDS-PAGE..
Source	Escherichia Coli
Expiration Date	12 months from date of receipt.